

Characteristics of Available Continuous Glucose Monitors

	Type	How Long Can Sensors Be Worn?	Are There Alarms?	Is Finger-Stick Testing Required for Calibration?
Medtronic MiniMed	Personal	Enlite sensors: 6 days	Yes	Yes
Medtronic iPro2 Professional	Diagnostic	Enlite sensors: 6 days	No	Yes
FreeStyle Libre	Personal	10 days	No	No
FreeStyle Libre Pro	Diagnostic	14 days	No	No
Dexcom G4 Platinum	Personal	7 days	Yes	Yes
Dexcom G5 Mobile	Personal	7 days	Yes	No
Dexcom G6	Personal	10 days	Yes	No